**Gebrauchsanweisung
Instructions for use
Instructions d'utilisation
Instrucciones de uso
Istruzioni per l'uso
Sposób użycia**

NADAL[®] Hb/Hp Complex (test cassettes)

REF 272011

Version 1.0

2010-11-04

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14. Referenzen

- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.
- and Mortality Weekly Report; (1995) 45 (5) : 107-110

1. Intended Use

Hb/Hp-Complex is a visual, immunochromatography dual quick double test for the qualitative detection of human hemoglobin in stool samples.

2. Clinical Significances

The test is intended for professional in vitro diagnostic use. The test was designed to supply evidence of bowel cancer and/or adenomas. Bowel cancer is one of the most frequently diagnosed types of cancer and one of the most frequent causes of death from cancer. Bowel cancer can be detected early by means of tests for hidden blood in the stools and the mortality rate can therefore be reduced. The Hb/Hp complex plays an important role in the retrieval of hemoglobin from lysed erythrocytes and is relatively stable against acid and proteolytic degradation. This means that the Hb/Hp complex can be detected even after longer passages through the bowel. This means that blood mixed with larger intestinal polyps and higher carcinomas can also be detected. The Hb/Hp complex displays a significant increase in sensitivity with regard to the recognition of colorectal adenomas and carcinomas. The combined evidence of hemoglobin (Hb) and hemoglobin-haptoglobin complex is therefore an innovative and extremely sensitive method for the early recognition of colon carcinomas and larger adenomas.

3. Principle of the Test

The immunological NADAL[®] Hb/Hp complex is based on a immunochromatographic method, in which both Hb and Hb/Hp complexes are specifically recognized through specific antibody reactions.

For Hb:

The hemoglobin contained in the stool sample reacts with specific monoclonal antibodies that are attached to gold particles. This complex is distributed on the membrane and reaches the test band (T) on which the anti-hemoglobin is present.

For Hb/Hp:

An antibody attaches to hemoglobin and another antibody attaches to the haptoglobin. The result is that specific Hb/Hp complexes are measured and there is no cross-reactivity with hemoglobin. In the case of a positive result, the molecules from the stool sample loaded with gold-marked antibodies attach to the test band (T) and are visible by means of a pink/red colouration. In the case of a negative result, there are no hemoglobin molecules that can attach to the test band (T) as complexes and therefore, there can also be no colouration of the test band (T). If the control strip (C) turns red/pink in colour, this shows that the sample has been taken correctly and has migrated correctly. The test is therefore valid.

4. Reagents and Materials Supplied

- 25 individually packed test cassettes
- 25 sample-taking tubes
- 25 sets of instructions on taking stool samples
- 25 stool containers
- 1 set of test instructions

5. Additional Required Materials

- 1 stopwatch

6. Storage & Stability

The test can be stored at +2° to +30°C in a sealed bag until the expiration date.

7. Warnings and Precautions

- Valid only for in vitro diagnostic use and only for professional use.
- Do not use the test after the expiration date.
- Do not mix test collection containers with different batches.
- Only open the test packaging if the test is actually going to be performed.
- All patient samples should be handled as if they could carry disease.
- The extraction solution contains a small amount of sodium azide.
- Patients should follow the instructions for obtaining a sample exactly and should not take any sample during menstruation or in case of bleeding hemorrhoids.

It is recommended that the stool container enclosed should be used, in order to avoid mixing with blood from urine samples or interfering components from the toilet water.

8. Specimen Collection and Preparation

Store the tube at room temperature before use. Collect the stool using the stool container enclosed. Hold the tube vertically and open the sealing cap. Remove the spiral rod. Place the spiral rod at three different locations in the stool sample. Place the spiral rod with the stool sample into the tube and seal well. Label the well-sealed sample-taker with the name of the patient. Store the tube in as cool a place as possible overnight and protected from light. Stool samples should be tested in the buffer solution within 1 week. Within this period of time, it is possible to store, transport and test the sample at room temperature without needing to take further measures. If the sample is taken at the patient's home, the sample should be stored where protected from light and in a cool place.

9. Procedure of the Test

The test cassette and the patient's sample should be brought to room temperature (+20°C to +30°C) before testing.

Take the required number of test cassettes from the foil packaging only immediately before performing the test. Mark the cassette with the name of the patient or with another form of identification.

Carefully shake the stool sample tube to ensure that the stool sample mixes properly with the extraction solution.



Unscrew the protective cap. Take a paper towel and break the seal of the collection container with the rotary motion. Hold the collection container upright and add 3 drops of the solution into the round sample window (S) of the cassette.

Read the results after 5 minutes. Strong positive results may be evaluated even sooner. Evaluate the result a maximum of 15 minutes after dropping the solution into the round sample window.

10. Interpretation of the Results



Test result is evaluated as "positive", if at least one of the two test lines appears.

Positive: 2 pink-coloured lines appear; one in the control region (C) and one in the test region (T). If you are testing strongly positive samples, the intensity of the control line may be reduced. It is not recommended to compare the intensity of the lines.

Negative: 2 pink-coloured lines appear; one in the control region (C) and one in the test region (T).

Invalid: If no red line appears in either region, this is a sign that the test is not functioning properly, or that the test materials are not correct. In this case, repeat the test with a new test cassette or contact the manufacturer for technical support.

11. Quality Control

The test contains a procedural control A coloured line that appears in the control region (C) shows that each test was performed correctly.

It is easy for the background to become slightly yellowish in colour during testing, depending on the colour of the stool sample. This is acceptable, as long as the evaluation of the test results are not adversely affected.

12. Limitations

The test is used for the qualitative determination of human hemoglobin and hemoglobin-haptoglobin complexes in stool samples.

A positive test occurs if human hemoglobin or hemoglobin-haptoglobin complexes are present in the sample. In addition to intestinal bleeding, blood in stools may also have other causes, such as hemorrhoids or the mixing of the stool sample with urine containing blood.

Not all intestinal bleeding are the result of benign or malignant polyps. The information that you obtain by performing this test should be examined in conjunction with other clinical test methods.

Because carcinomas and polyps can bleed to different extents and intermittently, in the case of immunological verification procedures, it is also advisable to check various stool samples.

Urine and excessive dilution of samples with water from the toilet bowl can lead to false test results. It is therefore recommended that the stool container be used.

Stool samples should not be taken during menstruation or 3 days before or after, in the case of bleeding caused by constipation, bleeding hemorrhoids, or in the case of medicine being administered rectally. This could lead to false positive results.

Blood-thinning medications, such as for example, acetylsalicylic acid or coumarin, as well as iron preparations, can lead to bleeding not caused by tumors.

No special diet is necessary before the performance of the test.

Unlike test methods that detect only free hemoglobin, the Hb/Hp complex very sensitive in showing bleeding higher up in the intestinal passage.

13. Performance Characteristics

Analytical sensitivity

The limits of detection of the tests are 25ng of free or complex hemoglobin per milliliter of buffer solution.

Prozone effect: The Hb/Hp complex also works reliably at extremely high hemoglobin values (>500,000ng/mL).

Specificity

The NADAL[®] Hb/Hp complex is specifically for human hemoglobin or the hemoglobin-haptoglobin complex and has no cross-reactivity with hemoglobin from cattle, pigs, rabbits, horses, or sheep at a concentration of up to 0.5 mg/mL in the extraction buffer. The Hb/Hp complex also displays no cross-reactivity with Bilirubin, vitamin C, and horseradish peroxidase.

14. References

- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.
- and Mortality Weekly Report; (1995) 45 (5) : 107-110

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14. Références

- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods; Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.
- and Mortality Weekly Report; (1995) 45 (5) : 107-110

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14.Referencias

Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.
and Mortality Weekly Report; (1995) 45 (5) : 107-110

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nal von minden GmbH