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**PRÄZITHERM®**  
**High temperature hotplate**

made of TITANIUM up to 600°C  
 for continuous operation


with PID-Temperature Controller  
 or  
 Programmer

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## High temperature hotplate PZ28-3T/D made of Titanium up to 600°C

This high temperature hotplate reaches its max. temperature of 600°C within about 20 minutes. The upper plate consists of TITANIUM, overall dimension 280 x 200 mm. The marked area of constant temperature is 230 x 160 mm. Due to a special insulation the device can be placed even on surfaces sensitive to heat. Metal sheets all around keep the outer casing temperature to a minimum.

The insulated hinged lid ensures almost no heat loss and quick heat-up time. Between the closed lid (inner surface CERAN®) and the heating plate remain 26 mm, enabling an intensive heating of plain objects at high precise temperatures. The whole casing is made of stainless steel.

The temperature regulator is connected by two plugged cables.

PZ28-3T(D) is available either with a temperature controller or with a programmer. The temperature controller keeps the device constant to an adjusted temperature within the whole range.

The programmer allows the control of certain programs including phases of heating up, constant temperature and cooling.

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### TEMP. CONTROLLER TR28-3T

This microprocessor-controlled temperature regulator is suitable for the high temperature hot plate. After adjustment of the temperature setpoint, the PID controller keeps it constant (control deviation 2 K or less).

### PROGRAMMER PR5-3T

The microprocessor-controlled ramp regulator offers a comprehensive time/temperature regulation for the high temp. hot plate.

Pressing the START key activates a program of up to 5 ramps. Max. 5 ramps (consisting of ramp time, constant temperature and time of constant temperature) can be programmed. Three program sets can be loaded into the non-volatile memory.

A programmed lead time enables the activation of the start up to 99 hours in advance.

Every temperature and time parameter can be selected with its own key. Consequently the setting is easy and clear. It is possible to pass over single ramps by setting the according parameter to "0". Every set value can be controlled during the running program. Data-Back-up is made by an EEPROM. The device is equipped with a change-over switch for continuous control. A buzzer signal sounds at the end of the program.

### Equipment of PR5-3T and TR28-3T

2-poled main switch; individual adjustment of proportional range, derivative/integrating factor and cycle time; control range limitation (to avoid the setting of values not useful or harmful to an application); contact maker input, actual value correction, junction bush for high temperature hot plate, wattage power control 1...100%, uninterrupted switching.

Both regulators show the actual value constantly. Furthermore, the programmer shows the remaining time of a ramp. When the SET key is pressed, the display changes to show the setpoint.

On request the devices are available with a temperature limiter to avoid overheating.

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### Further Fabrications

CERAN® hotplates up to 500°C - Precision hotplates - Precision hotplate made of cast iron up to 450°C - Stand hotplates with supporting rod - Rapid incinerator - Multipurpose extremely deep heating bath - Sandbaths - Series heating units - Multipurpose heating baths - Water- and Oil-baths - Development - Special orders.