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**BASELINE® EVALUATION
INSTRUMENTS**

HYDRAULISCHER HANDDYNAMOMETER

Fabrication
Enterprises
Incorporated

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Manufacturer and Master Distributor of
Physical Therapy and Rehabilitation Products

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
Virgil Mathiowetz, PhD, OTR
Associate Professor
Department of Occupational Therapy
College of St. Catherine
2004 Randolph Avenue
St. Paul, MN 55105-1794

March 18, 1993

Mr. Elliott Goldberg,
Marketing Director
Fabrication Enterprises Inc.
Trent Building
South Buckout Street
Irvington, NY 10533

Dear Mr. Goldberg,

Recently, I completed the study to determine whether the Baseline and Jamar hydraulic dynamometers can be used interchangeable. A draft of the report has been completed and sent to you. In the summary, I concluded that, "The data from this study suggest that the Jamar and Baseline hydraulic hand dynamometers measure equivalently for practical purposes. As a result, individuals using the Baseline Dynamometer are justified in using the normative data, which was collected with the Jamar dynamometer (Mathiowetz et al., 1985; 1986)." This conclusion assumes that the same standard procedures are followed as were used in the original normative data studies.

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

REFERENCES:

1. Gill D., Reddon J., Renney C., Stefanyk W.: Hand Dynamometer: Effects of Trials and Sessions. *Perpetual and Motor Skills* 61: 195-8, 1985.
2. Everett P., Silks F.: The Relationship of Grip Strength to Stature, Somatotype Components, and Anthropometric Measurements of the Hand. *The Research Quarterly* 23: 161-6, 1952.
3. Mathiowetz V., Federman S., Wiermer D.: Grip and Pinch Strength: Norms for 6 to 19 Year Olds. *The American Journal of Occupational Therapy* 40: 705-11, 1986.
4. Mathiowetz V., Donahoe L., Renells C.: Effect of Elbow Position on Grip and Key Pinch Strength. *The Journal of Hand Surgery* 10A: 694-7, 1985.
5. Mathiowetz V., Kashman N., Volland G., Weber K., Dove M., Rogers S.: Grip and Pinch Strength: Normative Data for Adults. *Archives of Physical Medicine and Rehabilitation* 66: 69-74, 1985.

* charts generated from data published in Mathiowetz's article "Grip and Pinch Strength: Normative Data for Adults", *Archives of Physical Medicine and Rehabilitation* 66: 69-74, 1985.