Art.Nr. 2192505
2192518
2192521
2192606
2192619
2192622
2192811
2192824
2192925
2193029
2193221

SISSEL[®] SITFIT[®] /PLUS

- ▮ Gebrauchsanweisung
- ▮ Instructions for use
- ▮ Mode d'emploi
- ▮ Gebruiksaanwijzing
- ▮ Instrucciones de uso
- ▮ Istruzioni per l'uso
- ▮ Bruksanvisning
- ▮ 사용설명서

SISSEL® SITFIT®

SISSEL® SITFIT® Plus (Ø 37 cm)

SISSEL® SITFIT®
(Ø 33 cm & Ø 36 cm)

SISSEL® SITFIT®

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SISSEL® SITFIT®

**THE SAFE AND COMFORTABLE WAY TO ACTIVELY IMPROVE YOUR SITTING POSTURE.
A CONVENIENT AND SIMPLE ALTERNATIVE TO A REGULAR GYM BALL, ERGONOMIC,
DYNAMIC AND PRACTICAL.**

The SISSEL® SITFIT® is a portable, air-filled sitting pad and therefore not just a dynamic sitting aid but also a training device for exercises while standing or sitting, in supine, lateral or prone position.

1 The SISSEL® SITFIT® is mainly a dynamic sitting aid, when placed under the feet the instable pad enables you to do leg exercises during long sitting periods.

2 Dynamic Sitting strengthens the stabilizing muscles of the back and therefore relieves the spine. Especially people who sit on the job use the SISSEL® SITFIT® with their office chair or at home. In addition it can be used to target restless sitting, for example at school.

The SISSEL® SITFIT® is skidproof and secure, which makes it a good choice for people who reject a gym ball as dynamic sitting device because of its size and possible instability.

Indications

- | Postural back pains (lumbar, thoracic, cervical)
- | Head-aches caused by postural deficiencies
- | Muscular disbalance of the torso
- | Prevention of back pains caused by prolonged sitting

Contra indications

- | Instable fractures of spine or pelvis
- | Open wounds
- | Any ailments for which the treating physician advises against the use of an air-filled pad

Sitting on the SISSEL® SITFIT® should initially be limited to only a few minutes to help the body getting used to it. This time can be gradually increased by a few minutes every day. Sitting should be altered between a regular chair and dynamic sitting on the SISSEL® SITFIT®.

Exercises to strengthen muscles of back, pelvis floor, legs and arms are possible while standing, sitting, in supine, lateral and prone position. The hardness of the SISSEL® SITFIT® can be individually adjusted with the help of a pump (scope of delivery with the SISSEL® SITFIT® Plus).

3 The needle pump permits the addition or removal of air whenever desired. Just screw the tip of the needle onto the pump and carefully insert into the valve. Do this in a careful manner in order to avoid any damage of the SISSEL® SITFIT®'s outer membrane.

Pressing your knee or lower arm onto the SISSEL® SITFIT® is helpful because it creates a bulge around the valve (this avoids a puncturing of the product). Now insert the pump and add air. To deflate simply the needle into the valve, hold slightly sideways until the desired amount of air has escaped.

4 Depending on the surrounding temperature the SISSEL® SITFIT® might feel a bit cool at first but it adjusts to your body temperature after a few minutes. Our attractive velour covers, available in many colors, create a pleasant feel instantly.

The SISSEL® SITFIT® is made in Germany. The PVC used for the product is free of phthalates. Instead we use CITROFOL® BII as plasticizer, a high-quality PVC softener based on citric acid. It is pharmacologically and eco-toxicologically harmless and is used by the food industry.

Note

- | Keep the application area free of pointy or sharp edges.
- | Avoid any difficult exercises or experiments with increased risk.
- | CITROFOL® BII as raw material features a high tolerance to light. However, to maintain the color brilliance we recommend not to expose your SISSEL® SITFIT® to direct sunlight for an extended time.
- | To clean wipe, your SISSEL® SITFIT® with a moist towel, water and using a mild household cleaner or standard disinfectant.
- | Weight tested up to 138 kg / approx. 300 lbs

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SISSEL® GmbH