Gebrauchsanweisung
Instructions for use
Instructions d'utilisation
Instrucciones de uso
Istruzioni per l'uso
Sposób użycia

nal von minden **Drug-Screen[®] Pipette test**

Version 1.0

2011-01-24

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7. Literatur

- Aniline O., Pittes, F. N., Phencyclidine (PCP): A review and perspectives. CRC Crit. Rev. Toxicol, 1982, 10, 145-177.
- Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
- Thomas L. eds., Labor und Diagnose, 6. ed., TH-Books Verlagsgesellschaft, Frankfurt, 2005
- Urine Testing for Drugs of Abuse, National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.
- Ellenhorn, M.J. and Barceloux, D.G. Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69,11970-11979, 1988.
- Gilman, A. G., and Goodman, L. S. The Pharmacological Basis of Therapeutics, eds. MacMillan Publishing, New York NY, 1980.
- Gorodetzky, C. W., Detection of Drugs of Abuse in Biological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring - Verlag, 1977.
- Greenblatt, D.J., Shader, R.I. Benzodiazepines in Clinical Practice. New York: Raven Press, 1974.
- Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.
- Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
- McBay, A. J., Clin. Chem. 33, 33B-40B, 1987.

Ref: 2011-01-24 SvM

1. Intended Use

The single-sandwich immunoassays are for the qualitative determination of various drugs and drug metabolites in human urine. The test is for professional use only (medical professionals).

2. Procedure of the Test

- 1) Remove the test from the packaging. Bring urine samples to room temperature before testing them.
- 2) Transfer 3 drops of urine (about 150 µl) into each sample opening of the test cassette. Make sure that the urine does not come into contact with other areas of the test.
- 3) Read the results after 5 minutes.


Transfer the 3 drops of urine here.

3. Interpretation of the Results


The test zone (T) will detect the presence of the specified drugs in the urine. The control zone is the area labeled (C)

> With a **negative result**, a drug-free urine sample, a pink stripe will appear in test zone (T) as well as a pink strip in the control area (C)

All of the parameters are negative.

> With a **positive result** a pink strip will only appear in the control area (C)

The parameter for cocaine (COC) is positive.

If no strip appears in the control zone (C) example test is invalid!

4. Note

The intensity of the colour stripe can vary as to what parameters are being tested. Only when no strip appears in the test zone than the test is positive. Positive results should be confirmed with another method (GC/MS).

Burprenorphine can be reliably detected from a daily dose of 2 mg.

5. Warnings and Precautions

- Only to be used for *in vitro* diagnostics.
- Only to be used by medical professionals.
- Do not use the test after the expiry date; the tests are no longer valid.
- Use necessary precautions when handling the urine samples
- Store the tests at 2-30°C, be sure not to freeze.
- Do not use the tests that have damaged packaging.
- Use test immediately after removal from protective packaging.
- Tests are good for only 1 use.

6. Cut-offs

Parameter	Drug / Metabolite**	Cut-off (ng/ml)
AMP (Amphetamine)*	Amphetamine	1000
AMP (Amphetamine)*	Amphetamine	500
AMP (Amphetamine)*	Amphetamine	300
BAR (Barbiturate)	Secobarbital	300
BZD (Benzodiazepine)*	Oxazepam	300
BZD (Benzodiazepine)*	Oxazepam	200
BZD (Benzodiazepine)*	Oxazepam	100
BUP (Buprenorphine)	Buprenorphine-β3-D-	10
COC (Cocain)*	Benzoylcegonine	300
COC (Cocain)*	Benzoylcegonine	100
COT (Cotinine)	Cotinine	200
EDDP (Methadone Metabolite)	2-Ethylidene-1,5-Dimethyl-3,3-Diphenylpyrrolidine	100
FYL (Fentanyl)	Fentanyl	10
KET (Ketamine)	Ketamine	1000
MDA (Amphetamine)	Methylenedioxyamphetamine	500
MDMA/ XTC (Ecstasy)*	3,4-Methylenedioxy-	1000
MDMA/ XTC (Ecstasy)*	3,4-Methylenedioxy-	500
MET (Methamphetamine)*	Methamphetamine	1000
MET (Methamphetamine)*	Methamphetamine	500
MET (Methamphetamine)*	Methamphetamine	300
MOR 2000	Morphine	2000
MOR 300 (Morphine/Heroin)*	Morphine	300
MOR 100 (Morphine/Heroin)*	Morphine	100
MQL (Methaqualone)	Methaqualone	300
MTD (Methadone)	Methadone	300
OXY (Oxycodone)	Oxycodone	100
PCP (Phencyclidine)	Phencyclidine	25
PCM (Paracetamol)	Paracetamol/Acetaminophen	5000
PPX (Propoxyphene)	Norpropoxyphene	300
TCA (Tricyclic Antidepress.)	Nortriptyline	1000
THC (Marihuana)*	11-nor-Δ9-THC-9-COOH	150
THC (Marihuana)*	11-nor-Δ9-THC-9-COOH	50
THC (Marihuana)*	11-nor-Δ9-THC-9-COOH	25
TML (Tramadol)	Tramadol	100

* The cut-off level can vary

7. Literature

- Aniline O., Pittes, F. N., Phencyclidine (PCP): A review and perspectives. CRC Crit. Rev. Toxicol. 1982, 10, 145-177.
- Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
- Thomas L. eds., Labor und Diagnose, 6. ed., TH-Books Verlagsgesellschaft, Frankfurt, 2005
- Urine Testing for Drugs of Abuse, National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.
- Ellenhorn, M.J. and Barceloux, D.G. Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.
- Gilman, A. G., and Goodman, L. S. The Pharmacological Basis of Therapeutics, eds. MacMillan Publishing, New York NY, 1980.
- Gorodetzky, C. W., Detection of Drugs of Abuse in Biological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring - Verlag, 1977.
- Greenblatt, D.J., Shader, R.I. Benzodiazepines in Clinical Practice. New York: Raven Press, 1974.
- Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.
- Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
- McBay, A. J., Clin. Chem. 33, 33B-40B, 1987.

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7. Riferimenti Bibliografici

- Aniline O., Pittes, F. N., Phencyclidine (PCP): A review and perspectives. CRC Crit. Rev. Toxicol. 1982, 10, 145-177.
- Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
- Thomas L. eds., Labor und Diagnose, 6. ed., TH-Books Verlagsgesellschaft, Frankfurt, 2005
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.
- Ellenhorn, M.J. and Barceloux, D.G. Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.
- Gilman, A. G., and Goodman, L. S. The Pharmacological Basis of Therapeutics, eds. MacMillan Publishing, New York NY, 1980.
- Gorodetzky, C. W., Detection of Drugs of Abuse in Biological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring - Verlag, 1977.
- Greenblatt, D.J., Shader, R.I. Benzodiazepines in Clinical Practice. New York: Raven Press, 1974.
- Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.
- Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
- McBay, A. J., Clin. Chem. 33, 33B-40B, 1987.

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7. Bibliografia

- Aniline O., Pittes, F. N., Phencyclidine (PCP): A review and perspectives. CRC Crit. Rev. Toxicol, 1982, 10, 145-177.
- Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
- Thomas L. eds., Labor und Diagnose, 6. ed., TH-Books Verlagsgesellschaft, Frankfurt, 2005
- Urine Testing for Drugs of Abuse, National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.
- Ellenhorn, M.J. and Barceloux, D.G. Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.
- Gilman, A. G., and Goodman, L. S. The Pharmacological Basis of Therapeutics, eds. MacMillan Publishing, New York NY, 1980.
- Gorodetzky, C. W., Detection of Drugs of Abuse in Biological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring - Verlag, 1977.
- Greenblatt, D.J., Shader, R.I. Benzodiazepines in Clinical Practice. New York: Raven Press, 1974.
- Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.
- Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
- McBay, A. J., Clin. Chem. 33, 33B-40B, 1987.

1. Detailed list of cut-offs

The following structurally related substances resulting from the concentration following a positive result.

Amphetamine

Amphetamine -Test with Cut-off 300 ng/ml (AMP300).

Drug	Concentration (ng/ml)
D-Amphetamine	300
3,4-Methylenedioxy-amphetamine (MDA)	500
3,4-Methylenedioxy-methamphetamine (MDMA)	>20.000
L-Amphetamine	>50.000

Amphetamine -Test with Cut-off 500 ng/ml (AMP500).

Drug	Concentration (ng/ml)
D-Amphetamine	500
L-Amphetamine	15.000
(+/-)3,4-methylenedioxyamine (MDA)	700
MDMA	12.000

Amphetamine-Test with Cut-off 1000 ng/ml (AMP1000)

Drug	Concentration (ng/ml)
D-Amphetamine	1.000
L-Amphetamine	>50.000
D-Methamphetamine	>20.000
L-Methamphetamine	>20.000
3,4-Methylenedioxyamphetamine (MDA)	2.400
3,4-Methylenedioxy-methamphetamine (MDMA)	>20.000
3,4-Methylenedioxyethylamphetamine (MDEA)	>100.000
Paramethoxyamphetamine (PMA)	1.000

Barbiturate

Barbiturat-Test with Cut-off 300 ng/ml (BAR)

Drug	Concentration (ng/ml)
Secobarbital	300
Phenobarbital	300
Butalbital	3.000
Allobarbital	5.000
Alphenal	625
Amobarbital	600
Aprobarbital	600
Hexobarbital	>100.000
Butobarbital	75
Pentobarbital	300

Buprenorphine

Buprenorphine -Test with Cut-off 10 ng/ml (BUP).

Drug	Concentration (ng/ml)
Buprenorphine	10
Buprenorphine-3- β -d-glucuronide	10
Nor-Buprenorphine	>1.000
Nor-Buprenorphine-3- β -d-glucuronide	>1.000

Benzodiazepine

Benzodiazepin-Test with Cut-off 100 ng/ml (BZD100)

Drug	Concentration (ng/ml)
Oxazepam	100
a-Hydroxyalprazolam	625
Alprazolam	63
Bromazepam	250
Chlordiazepoxid	1.250
Clobazam	31
Clonazepam	625
Clorazepat	1.670
Delorazepam	625
Desalkflurazepam	125
Diazepam	125
Estazolam	1.250
Flunitrazepam	125
Lorazepam	313
Lormetazepam	750
Midazolam	2.500
Nitrazepam	12.500
Norchlordiazepoxid	63
Nordiazepam	250
Temazepam	63
Triazolam	750

Benzodiazepine-Test with Cut-off 200 ng/ml (BZD200)

Drug	Concentration (ng/ml)
Oxazepam	200
a-Hydroxyalprazolam	1.250
Alprazolam	188

Bromazepam	500
Chlordiazepoxide	2.500
Clobazam	63
Clonazepam	1.250
Clorazepat	3.330
Delorazepam	1.250
Desalkflurazepam	400
Diazepam	250
Estazolam	2.500
Flunitrazepam	250
Lorazepam	1.250
Lormetazepam	1.250
Midazolam	7.500
Nitrazepam	25.000
Norchlordiazepoxid	125
Nordiazepam	500
Temazepam	125
Triazolam	1.250

Benzodiazepine-Test with Cut-off 300 ng/ml (BZD)

Drug	Concentration (ng/ml)
Oxazepam	300
a-Hydroxyalprazolam	2.500
Alprazolam	250
Bromazepam	625
Chlordiazepoxid	3.750
Clobazam	100
Clonazepam	2.500
Clorazepat	5.000
Delorazepam	2.500
Desalkflurazepam	500
Diazepam	300
Estazolam	5.000
Flunitrazepam	375
Lorazepam	1.500
Lormetazepam	2.000
Midazolam	10.000
Nitrazepam	37.500
Norchlordiazepoxid	250
Nordiazepam	750
Temazepam	150
Triazolam	1.875

Cocaine

Cocaine-Test with Cut-off 100 ng/ml (COC100)

Drug	Concentration (ng/ml)
Cocaine	3.000
Benzoyllecgonine	100
Ecgonine-Methylester	>100.000
Ecgonine	>40.000

Cocaine-Test mit Cut-off 300 ng/ml (COC300)

Drug	Concentration (ng/ml)
Cocaine	1.000
Benzoyllecgonine	300
Ecgonine-Methylester	>100.000
Ecgonine	>40.000

Cotinine-Test

Cotinine-Test mit Cut-off 200 ng/ml (COT)

Drug	Concentration (ng/ml)
Cotinine	200
(-)-Nicotine	>50.000

EDDP/Methadone-Metabolites

EDDP-Test with Cut-off 100 ng/ml (EDDP)

Drug	Concentration (ng/ml)
2-ethyliden-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	100
Doxylamine	>1.000.000
Methadone	>10.000
Methadol	>10.000

Fentanyl-Metabolites

Fentanyl-Test with Cut-off 10 ng/ml (Fentanyl)

Drug	Concentration (ng/ml)
Fentanyl-Metabolite	10

Ketamine

Ketamine-Test with Cut-off 1000 ng/ml (Ketamine)

Drug	Concentration (ng/ml)
Ketamine	1.000
Norketamine	2.000
Phencyclidine	>10.000
Dextromethorphan	1.000
Lidocaine	>100.000
Promethazine	20.000

d-Amphetamine	>100.000
Methadone	25.000
Methadol	50.000
Pethidine	12.500
Methamphetamine	12.500
MDMA	>100.000
MDA	
MDA-Test with Cut-off 500 ng/ml (MDA)	
Drug	Concentration (ng/ml)
(+/-)3,4-Methylenedioxyamphetamine (MDA)	500
3,4-Methylenedioxyethylamphetamine	300
MDMA/XTC	
MDMA-Test with Cut-off 500ng/ml (MDMA/XTC500)	
Drug	Concentration(ng/ml)
3,4-Methylendioxy-methamphetamine (MDMA)	500
3,4-Methylendioxyamphetamine (MDA)	1.000
3,4-Methylendioxyethylamphetamine (MDEA)	300
D-Amphetamine	>100.000
D-Methamphetamine	>100.000
Paramethoxyamphetamine (PMA)	5.000
MDMA/XTC-Test with Cut-off 1000 ng/ml (MDMA/XTC1000)	
Drug	Concentration(ng/ml)
3,4-Methylendioxy-methamphetamine (MDMA)	1.000
3,4-Methylendioxyamphetamine (MDA)	2.000
3,4-Methylendioxyethylamphetamine (MDEA)	600
D-Amphetamine	>100.000
D-Methamphetamine	100.000
L-Methamphetamine	>100.000
Methamphetamine	
Methamphetamine-Test with Cut-off 300 ng/ml (MET300).	
Drug	Concentration (ng/ml)
(+)-Methamphetamine	300
D-Amphetamine	>40.000
Chloroquine	5.000
(+/-)-Ephedrine	>100.000
L(-)-Methamphetamine	10.000
Mephentermine	50.000
(+/-)3,4-Methylenedioxy-methamphetamine (MDMA)	1.000
(+/-)3,4-Methylenedioxy-ethylamphetamin e(MDEA)	10.000
Procaine	35.000
β-Phenylethylamine	35.000
Ranitidine	35.000
Methamphetamine-Test with Cut-off 500 ng/ml (MET500).	
Drug	Concentration (ng/ml)
(+)-Methamphetamine	500
D-Amphetamine	>50.000
Chloroquine	>50.000
(+/-)-Ephedrine	>50.000
L(-)-Methamphetamine	>50.000
Mephentermin	>50.000
(+/-)3,4-Methylenedioxy-methamphetamine (MDMA)	700
Procaine	20.000
□-Phenylethylamine	30.000
Ranitidin	>50.000
Methamphetamine-Test with Cut-off 1000 ng/ml (MET1000)	
Drug	Concentration (ng/ml)
D-Methamphetamine	1.000
D-Amphetamine	>40.000
Chloroquine	10.000
(+/-)-Ephedrine	>100.000
L-Methamphetamine	15.000
3,4-Methylendioxy-amphetamine (MDA)	>10.000
3,4-Methylendioxy-methamphetamine (MDMA)	2.000
3,4-Methylendioxy-ethylamphetamine (MDEA)	20.000
Procaine	100.000
Methaqualone	
Methaqualone-Test with Cut-off 300 ng/ml (MQL)	
Drug	Concentration (ng/ml)
Methaqualone	300

Mecloqualone	500
3'-Hydroxy methaqualone	500
4'-Hydroxy methaqualone	500
2'-Hydroxy methaqualone	3.000
Amitriptyline	50.000
Carbamazepine	20.000
Nortriptyline	50.000
Phenytoine	40.000
Primidone	24.000
Theophylline	40.000
Methadone & Derivats	
Methadone-Test with Cut-off 300 ng/ml (MTD)	
Drug	Concentration (ng/ml)
Methadone	300
Methadole	1.000
2-Ethyliden-1,5-Dimethyl-3,3-Diphenylpyrolidine (EDDP)	>40.000
Doxylamine	>40.000
Opiate/Morphine	
Morphine-Test with Cut-off 100 ng/ml (MOR100).	
Drug	Concentration (ng/ml)
Morphine	100
Codeine	100
Diacetylmorphine	100
Ethylmorphine	100
Hydromorphone	500
Hydrocodone	500
Meperidine	>100.000
6-Monoacetylmorphine	100
Morphin-3-β-d-glucuronide	2.000
Oxycodone	>20.000
Oxymorphone	>20.000
Promethazine	>250.000
Rifampicin	8.400
Thebain	8.400
Trimipramine	>20.000
Morphine-Test with Cut-off 300 ng/ml (MOR300)	
Drug	Concentration (ng/ml)
Morphine	300
Codeine	300
Diacetylmorphine (Heroine)	300
Ethylmorphine	300
Hydromorphone	1.500
Hydrocodone	1.500
Merperidine	>100.000
6-Monoacetylmorphine	300
Morphin-3- β-d-glucuronide	6.000
Oxycodon	>20.000
Oxymorphon	>20.000
Promethazin	>250.000
Rifampicin	25.000
Thebain	2.50
Trimipramine	>20.000
Morphine-Test with Cut-off 2000 ng/ml (MOR2000).	
Drug	Concentration (ng/ml)
Morphine	2.000
Codeine	2.000
Diacetylmorphine (Heroine)	2.000
Ethylmorphine	600
Hydromorphone	15.000
Hydrocodone	15.000
Merperidine	>100.000
6-Monoacetylmorphine	5.000
Morphin-3-β-d-glucuronide	10.000
Oxycodon	>20.000
Oxymorphone	>20.000
Rifampicin	>50.000
Thebain	20.000
Morphine	2.000
Oxycodone	
Oxycodone-Test with Cut-off 100 ng/ml (OXY).	
Drug	Concentration (ng/ml)
Oxycodone	100
Oxymorphone	500
Codeine	>100.000
Diacetylmorphine (Heroine)	>100.000
Ethylmorphine	10.000
Hydrocodone	10.000
Hydromorphone	10.000
Merperidine	>100.000
6-Monoacetylmorphine	>100.000
Morphin-3-beta-D-glucuronide	>250.000
Thebain	>100.000

Phencyclidine	
PCP-Test with Cut-off 25 ng/ml (PCP)	
Drug	Concentration(ng/ml)
Phencyclidine	25
Tenocyclidine	2.000

Paracetamol	
PCM-Test with Cut-off 5000 ng/ml (PCM)	
Drug	Concentration (ng/ml)
Paracetamol	5.000

Propoxyphene	
PPX-Test with Cut-off 300 ng/ml (PPX)	
Drug	Concentration (ng/ml)
D-Propoxyphene	300
D-Norpropoxyphene	300

TCA	
TCA-Test with Cut-off 1000 ng/ml (TCA)	
Drug	Concentration(ng/ml)
Nortriptyline	1.000
Amitriptyline	1.000
Chlorpromazine	3.500
Clomipramine	10.000
Cyclobenzaprine	1.500
Desipramine	500
Diphenylraminee	20.000
Doxepine	1.000
Imipramine	800
Nordoxepine	1.000
Opi Pramol	4.000
Protriptylin	3.000
Doxepin	1.000
Perphenazine	25.000
Promazin	200
Promethazin	40.000
Protryptilin	3.000
Trimipramine	2.500

THC & Derivats	
THC-Test with Cut-off 25 ng/ml (THC25).	
Drug	Concentration (ng/ml)
11-nor- Δ^9 -THC-9-COOH	25
11-nor- Δ^8 -THC-9-COOH	25
11-hydroxy- Δ^9 -Tetrahydrocannabinol	>100.000
Δ^8 -Tetrahydrocannabinol	10.000
Δ^9 -Tetrahydrocannabinol	10.000
Cannabinol	15.000
Cannabidiol	>100.000

THC-Test with Cut-off 50 ng/ml (THC50)	
Drug	Concentration(ng/ml)
11-nor Δ^8 -THC-9-COOH	50
11-nor- Δ^9 -THC-9-COOH	50
11-hydroxy- Δ^9 -Tetrahydrocannabinol	>100.000
Δ^8 - Tetrahydrocannabinol	15.000
Δ^9 - Tetrahydrocannabinol	15.000
Cannabinol	20.000
Cannabidiol	>100.000



THC-Test with Cut-off 150 ng/ml (THC150)	
Drug	Concentration (ng/ml)
11-nor- Δ^8 -THC-9-COOH	150
11-nor- Δ^9 -THC-9-COOH	150
11-hydroxy- Δ^9 -Tetrahydrocannabinol	>100.000
Δ^8 -Tetrahydrocannabinol	80.000
Δ^9 -Tetrahydrocannabinol	80.000
Cannabinol	>20.000
Cannabidiol	>100.000

Tramadol	
Tramadol-Test with Cut-off 100 ng/ml (Tramadol)	
Drug	Concentration (ng/ml)
Tramadol	100
N-desmethyl-tramadol	250
O-desmethyl-tramadol	10.000

With exception of the respective parameters positive reacting drugs and drug metabolites, which are listed above, all below listed compounds react negatively up to a concentration of 100 μ g/ml.

Acetaminophene	Guaiacol Glyceryl Ether
Acetone	Hämoglobin
Albumin	Ibuprofen
Ampicillin	(+/-)-Isoproterenol
Aspartam	Koffein
Aspirin	Kreatin
Atropin	Lidocain
Benzocaine	L-Phenylephrin
Beta-Phenylethylamine	N-Methyl-Ephedrin
Chinidine	(+)-Naproxen ([S]-6-Methoxy--Methyl-2-Naphthalenessigsäure)
Chloroquine	Oxalsäure
Chlorpheniramine	Penicillin-G (Benzylpenicillin)
4-Dimethylaminoantipyrine	Pheniramin
Dopamin (3-Hydroxytyramin)	Phenothiazin
(-)-Ephedrine	Procain
(+/-)-Ephedrine	Ranitidin
Ethanol	Sulindac
Furosemid	Tyramin
	Vitamin C (Ascorbinsäure)

Symbol	English	Deutsch	Français	Nederlands	Español	Italiano
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Gebruiksaanwijzing	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	European Conformity	CE-Konformitätskennzeichnung	Conformité aux normes européennes	Conform Europese richtlijnen	Conformidad europea	Conformità europea
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	In vitro diagnostisch gebruik	Para uso Diagnóstico in vitro	Per uso Diagnostica in vitro
	Catalogue number	Katalog-Nr.	Numéro de catalogue	Catalogus nummer	Número de catálogo	Numero di Catalogo
	Lot. No. / Batch code	Chargen-Nr.	Numéro de lot	Lot nummer	Número de lote	Numero di lotto
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Geschikt voor <n> tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungs-temperatur	Température de conservation	Opslag temperatuur	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Houdbaarheids datum	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabrikant	Fabricante	Fabbricante
Distributed by	Distributor	Vertreiber	Distributeur	Distributeur	Distribuidor	Distributore
	Content	Inhalt	Conditionnement	Inhoud	Contenido	Contenuto

Symbol	Polski	Suomi	Portugues	Dansk	Svenska	Ελληνικά
	Patrz: ulotka informacyjna	Katso käyttöohjeet	Consulte as instruções de utilização	Se brugsanvisning	Se bruksanvisningen	Εγχειρίδιο χρήστη
	Znak zgodności CE	CE-merkitty	Conformidade com as normas europeias	Europeaisk overensstemmelse	Europeisk överensstämmelse	Ευρωπαϊκή Συμμόρφωση
	Tylko do użytku in-vitro	In vitro-diagnostiikka	Diagnóstico in vitro	In vitro diagnostik	In vitro diagnostik	in vitro διαγνωστικό
	Numer katalogowy	Luettelonumero	Catálogo n.º	Katalognummer	Katalognummer	Αριθμός καταλόγου
	Numer serii	Eränumero	No do lote	Lot nummer	Batchnummer	Αριθμός Παρτίδος
	Wystarszajonce na "n" powtórzeń	Sisältää tarvikkeet "n" testiin		Indeholder tilstrækkeligt til "n" test	Innehåller tillräckligt till "n" tester	Περιεχόμενο επαρκές για «n» εξετάσεις
	Temperatura przechowywania	Säilytyslämpötila	Temperatura de conservação	Opbevarings-temperatur	Förvaringstemperatur	Θερμοκρασία αποθήκευσης
	Data ważności	Viimeinen käyttöpäivä	Prazo de validade	Udløbsdato	Utgångsdatum	Ημερομηνία λήξης
	Producent	Valmistaja	Fabricante	Producent	Tillverkare	Κατασκευαστής
Distributed by	Dystrybutor	Tukkumyyjä			Distributör	
	Zawartość	Sisältö	Conteúdo	Indhold	Innehåll	Περιεχόμενο

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nal von minden GmbH