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User manual Senso® Balance set of 2 16 cm and 20 cm diameter

Congratulations!

You are now the owner of a Senso® Balance set of 2, a product “made in Germany”.

This sports tool was produced in compliance with valid EU norms, it contains no heavy metals and is made of 100% recyclable plastic.

This tool is completing the Senso® Balance product family. A new feature is the improved contact with the ground. The air filling and the special Senso® knobs ensure an ideal stimulation of the proprioceptive senses and are especially suitable for all stepping and walking exercises. Ideal for perceptual training of children as well as adults.

Warning

- Avoid storage in direct sunlight or near heat sources, because otherwise the Senso Balance set of 2 might deform
- Avoid contact with sharp and pointed objects (e.g. broken glass, pointed rocks in shoe soles) because they can damage the Senso Balance set of 2. Do not throw the set!
- Use sport shoes, regular shoes or be barefoot while training on the Senso Balance set
- Prior to any use, the Senso Balance set should be checked for damage because otherwise a safe handling is not guaranteed
- The Senso Balance set has to rest on a horizontal, stable and dry ground, to make sure, that the set cannot slide away during the exercises.

Equipment details:

Each piece of the SensoBalance set is made of high-grade Rulon with a needle valve.

Size: 16 cm Ø, ca. 20 cm Ø

Load: 150 kg

Care:

Use mild standard detergents and disinfectants. Dissolvers as well as acidic or corrosive cleansers can harm the material. The SensoBalance set is waterproof.

How to inflate:

The air filling of the SensoBalance Set can be regulated by using a standard needle pump (e.g. ball pump) Therefore the Senso Balance Set is equipped with a needle valve.

Please perform changes to the air pressure inside the Senso Balance Set (inflation or bleeding) only, when the Senso Balance has room temperature and add at the most as much air, as is necessary to make the surface level. Always moisten the needle before carefully inserting it into the valve.