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Introduction

People of all ages have been sitting, exercising and relaxing on the large, round, colourful balls for half a century. Initially limited to medical applications, the large ball moved into nearly all private households thanks to preventive programmes initiated by health insurance providers and remains popular in sports and fitness, therapy and health, recreation and play to this day. Few exercise devices are this versatile: The ball promotes balance and coordination, proprioception and posture, relaxation, flexibility and strength.

In addition to its outstanding flexibility, the ball convinces young and old with its incredibly motivating nature – it makes every use a "rounded" experience. Try it for yourself!

Areas of application

In addition to therapeutic applications, the fitness balls are increasingly being used in sports and fitness courses such as "Fit-Ball-Aerobics" or "Core Stabilisation". As a training and seating aid, the balls challenge and promote the entire locomotor system as well as the sensory system. The unstable sitting surface trains conditional and coordinative abilities, promotes versatile balance reactions and thereby enhances day-to-day safety.

In therapy, ARTZT vitality® Fitness Balls are used for mobility, strength and coordination training in all aspects of orthopaedic and neurological rehabilitation. They are used in case of Parkinson's, central paralysis, stroke, head and brain injuries, carpal tunnel syndrome, rotator cuff injuries and tennis elbow. They are suitable for people of all ages – from foot and senior care to aerobics and sports. In short: The large balls are much more versatile as training devices than one would think at first glance.

Materials and characteristics

ARTZT vitality® Fitness Balls are made in Germany, can be recycled 10-15 times and do not contain any heavy metals hazardous to health. The balls are available in two versions:

- » Standard (NOT burst proofed)
- » Professional with ABS® = Anti-Burst-System® (load capacity up to 500 kg, burst proof up to 90 kg)

Thanks to the patented Crylon material blend (a leathery, extremely abrasion-resistant surface), the professional balls with ABS® offer a higher standard of safety than the standard balls. The ABS® ball cannot burst after punctiform damage but merely loses air slowly. This makes it especially suitable for rehabilitation applications or when training certain target groups (pregnant women, hip patients etc.). The balls are colour-coded in the following sizes:

Body height	Ø ball	
< 155 cm	45 cm	
156 - 165 cm	55 cm	
166 - 179 cm	65 cm	
> 179 cm	75 cm	

Handling

Not too soft and not too firm, the robust material assures a high level of seating comfort that is influenced to a large degree by the following three factors:

- » Inflation technique
- » Ball size
- » Perfect sitting position of the user

Inflation technique

(please be sure to also read the product insert!)

All ARTZT vitality® Fitness Balls can be inflated with a hand pump or much more quickly and easily with a compressor. A conventional bicycle pump is not suitable. In contrast to the standard ball, the professional ball requires a special process:

1. Inflate the professional fitness ball in stages until it reaches 80 to 90 percent of the maximum size. Then allow it to rest for at least two hours. Do not inflate at temperatures below 20°C.

2. Inflate the professional fitness ball to the required or maximum size. Since the leathery material needs time to stretch and expand, it may tear when subjected to improper treatment.

IMPORTANT: The ABS® system is activated after approximately 24 hours. For your safety, ensure you allow the professional ball to rest for 24 hours after inflation before using it for the first time. Please also observe this rest period after re-inflating the ball when necessary.

Just like a car tyre, the fitness balls also lose air after being stored for approximately 8-12 months. In order to restore the original condition, the balls must be re-inflated. During inflation and re-inflation, it is important to ensure that the balls are not overstretched – i.e. inflated beyond the maximum intended size. The simplest way to check the size is to use a piece of string. In order to do so, first determine the required length of string (diameter x 3.14) and then wrap the string around the ball. As soon as the piece of string tightens during inflation, the desired size has been reached.

Ball size

The size of the ball must be selected according to the body size. The following values are offered as a guideline:

Body height	Ø ball	
< 155 cm	45 cm	
156 - 165 cm	55 cm	
166 - 179 cm	65 cm	
> 179 cm	75 cm	

Perfect sitting position of the user

A “sitting test” should be conducted in order to select the right ball with confidence. In order to do so, the person sits on the middle of the ball. The hip / seat angle should be between 90 and 110 degrees in the sitting position. The thighs slope downwards slightly and the soles of both feet are fully on the floor.

In addition to the factors already mentioned above, the durability of and safe training with the fitness balls mainly depends on proper use:

- » The balls should always be checked for damage before active use.
 - » Remove sharp or pointed objects from the immediate vicinity of the training area.
 - » If balls get dirty after intensive use by various people, they can be cleaned easily using a damp cotton cloth. Soap and water may be used for cleaning very dirty balls. Only use the products after they have dried completely.
 - » Do not store the balls in direct sunlight since this affects the colour and material surface over the long term. Storage and use at extreme temperatures should also be avoided.
 - » Improper or careless use of the balls may lead to injuries. As a result, only safe and pain-free exercises should be completed.
- » The balls are not intended as toys for children! Children must be supervised during use!
 - » Please observe the amount of space required to train with a ball: Approximately 4 m² per participant.

Fitness-oriented preventive training focuses on improved muscular endurance, coordination, muscle development and the resulting effects on the body. Beginners should select 6-8 of the following exercises in order to achieve this goal. Each exercise should be repeated 15-20 times. The entire training sequence may be repeated 1-3 times. If the focus of training is more on endurance, the exercises may also be repeated up to 30 times. For personal training with a performance-oriented client, working with a low repeat count (6-12) is also possible. You should never experience EXTREMELY heavy strain while training.

Please observe the following principles:

- » Maintain correct posture and conduct movements in a slow, controlled manner in the initial position, during the exercise, and at the end. Breathing evenly is important.
- » Training should be pain-free. Take a break if an exercise becomes painful. Consult your doctor if pain is felt repeatedly.
- » Training with children requires constant supervision and guidance.
- » Study each exercise in detail and visualise the process.
- » ALWAYS exercise both the right and left side of the body even when the illustration only shows one side.
- » Always complete balance and equilibrium exercises at the beginning of your training session since these exercises require the most concentration.
- » If your focus of training is muscle development, a combined balance and strength training (with additional weights) on the ball is not recommended since you can only use 30-40% of your potential strength.

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The following exercises are used for warm-ups and to become accustomed to the new training device. Start slowly. Complete all 4 warm-up exercises for 1 to 2 minutes each. Then begin the main training session.

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Sit on the middle of the ball. The soles of both feet are flat on the floor.

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Now move your pelvis alternately to the right and to the left.

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From the same starting position, tilt your pelvis forward (rounding your back) and then back (slightly hollow back).

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Sit upright on the ball. Find a stable sitting position. Now alternately tap your
left and right heels.

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Sit upright on the ball. Find a stable sitting position. Now alternately tap your left and right toes.

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Once this comes easily, you can also put both feet on tiptoe at the same time.

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Now alternately raise the right and left leg. You can complete this exercise with your leg extended and / or flexed

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Shift your weight ahead slightly and stand up slowly.

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The knees remain slightly bent. Hold this position briefly and then slowly sit back down on the ball. (Caution: Always check to make sure the ball is still in the same place.)

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Stand behind the ball. The fingertips are touching the top of the ball lightly.

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Now slowly bend your knees (lower your buttocks as though you were sitting down) and simultaneously roll the ball ahead slightly as you extend your arms. Make sure your knees stay behind the tips of your feet.

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Kneel down beside the ball and support your upper body on the ball. Extend the top leg and lift it off the floor by 10-20 cm.

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The foot should not touch the floor when you lower your leg between repeats.

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Lie on your back and hold the ball between your knees. Now squeeze the ball firmly between your legs. Hold briefly and then release.

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Kneel in front of the ball. Your hands are resting on the front of the ball.

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Now push off the floor with your feet and roll over the ball into a push-up position. Hold this position for a few seconds and then return to a kneeling position.

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Now push off the floor with your feet and roll over the ball into a push-up position. Do a push-up.

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"Pull" the ball closer to your hands by assuming a kneeling position (on the ball).

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From there, return to the push-up position and then to the kneeling position on the floor.

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Kneel in front of the ball and support your upper body on it. The hands are held loosely behind the head.

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Lift your upper body. Keep your gaze directed ahead and down.

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Once this exercise comes easily, you can also rotate your upper body alternately to the right and left as you lift it.

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Kneel in front of the ball and support your upper body on it. The arms are extended towards the floor in front of the ball.

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a) Slightly raise your upper body while moving your extended arms to the sides of your body (tighten the shoulder blades). Keep your gaze directed ahead and down.

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b) Slightly raise your upper body while moving your extended arms up. Keep your gaze directed ahead and down.

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Lie on the exercise ball on your stomach. Your hands and feet lightly touch the floor.

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Now extend the right leg and the left arm simultaneously. Hold this position for a few seconds and then change sides.

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Kneel down beside the ball and support your upper body on the ball. The left arm is bent, the hand rests lightly on the temple, the upper leg is extended.

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Slightly raise the upper body (elbow pulls towards the foot). Hold this position briefly and then lower your upper body.

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Sit down on the ball. The hands are resting lightly beside the thighs.

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Now slowly walk your feet forward while rolling the upper body over the ball. Briefly hold this position and then return to a sitting position.

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Lie on your back as described for the rollover and place both hands lightly behind your head.

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Now raise the upper part of your back off the ball (slightly lift your upper body). Your gaze is directed towards the ceiling and the elbows point out. Hold this position briefly and then lower your upper body. Once this exercise comes easily, you can also rotate your upper body alternately to the right and left as you lift it.

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Lie on your back as described for the rollover. The arms are extended and point towards the ceiling.

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Now alternately point towards the ceiling with the fingertips of the left and right hands. The upper body and / or the shoulder blades lift off the ball slightly at the same time.

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Dorsal position. The legs are bent. Both feet are on the floor. The ball is held in front of the upper body with the arms extended.

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Now slowly roll up into a sitting position by pushing the ball forward and up. In the sitting position (not fully upright, but with the back slightly bent), throw the ball up and catch it; then return to the initial position.

Ganzkörperkräftigung, Koordination | Full-body strength, coordination
Brücke | Bridge

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Dorsal position: The lower legs are supported on the ball,
the arms are extended alongside the body.

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Raise the hips so that the front of your body forms a straight line. Once you have found a stable position, you can extend one leg towards the ceiling (alternate legs).

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Then return to the position with the legs extended and finally to the initial position.

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Place both forearms on the ball and extend the legs by supporting yourself on the ball. The back of your body forms a straight line. Keep your balance in this position.

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Stretching

Regular stretching – especially after strength training – relaxes your musculature, shortens your regeneration time and reduces the risk of injury. Stretching may be done by holding a stretch or repeating it. When holding a stretch, you assume the position slowly until you feel the stretch (pulling but no pain) and then hold this position for 15-30 seconds. Repeat this type of stretch 2-3 times. With repeated stretching, also assume the stretch position slowly but then move into this stretch position slowly and in a controlled manner another 10-20 times (e.g. while exhaling). Both methods should be completed in a controlled manner (no sudden movements) and without causing pain. The following stretching exercises target the lower extremities, the back and the chest and shoulder musculature. To round out your stretching session, you should also complete a classic stretching exercise for your neck.

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Sit down to the left of the ball and extend your left leg backwards. Now slowly move the ball forward with the left leg (bounce). Switch legs.

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Sit on the middle of the ball. Extend the left leg and put your heel on the floor. Now bend the upper body forward slightly. Switch legs.

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Sit on the front section of the ball. Extend the right leg to the side and roll your body slightly to the left. Switch legs.

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Kneel down beside the ball and support your upper body on the ball. Extend your legs and raise your left arm straight above your head. The hand pulls slightly towards the floor. Change sides.

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Kneel approximately 50 cm in front of the ball. Place the palms of your hands on the ball and let your upper body “fall”. Now you can also move the ball slightly to the right and left using your arms.

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The children sit on the balls and hold them with their knees. Now they hop through the room forwards, backwards and sideways. They can also have a hopping race from one side of the room to the other.

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The children stand back to back and transport a ball through the room in this position. What directions can you move in without losing the ball? Can you also squat down?

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
The children roll over the ball into a push-up position (the stomach / hips are supported on the ball). Now they perform additional tasks such as clapping, swimming, moving to the right and left etc. with the hands / arms while in this position.

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The children lie on the ball so they are supported on their stomach. The feet touch the floor while the upper body and arms are raised (takeoff). Then the hands are supported on the floor while the legs are extended backwards and upwards (landing).

UPGRADE YOUR VITALITY!

ARTZT
vitality

Produktprogramm ARTZT vitality® | Product range ARTZT vitality® Medical Fitness - Made in Germany

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- **Flowin® Personal Training**
- **Gymstick™ Health**
- **BOSU® Balance Training**

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH. Unauthorized use is strictly prohibited.
© 2009 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

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© 2009 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.