








## HEINE XHL® Xenon Halogen Bulbs

	Instrument	2.5 V	3.5 V
	mini Fibrilux, mini 2000 Combi Lamp, mini 2000 Cliplamp, mini 1000 Cliplamp	<b>X-01.88.041</b>	
	alpha Ophthalmoscope, mini 2000 Ophthalmoscope, alpha Focalux, mini Miroflex, mini 2000 Focalux	<b>X-01.88.042</b>	
	alpha+ Ophthalmoscope , alpha+ Focalux	<b>X-01.88.042</b>	<b>X-02.88.101</b>
	alpha+ Ophthalmic Examination Light, alpha ophth. Examination Light – 2.5 V	<b>X-01.88.043</b>	<b>X-02.88.102</b>
	minilux Otoscope, mini 2000 Otoscope	<b>X-01.88.056</b>	
	Straight laryngeal mirror	<b>X-01.88.057</b>	<b>X-02.88.063</b>
	HSR 2 Retinoscope	<b>X-01.88.058</b>	<b>X-02.88.048</b>
	Conventional laryngoscope blades	<b>X-01.88.059</b>	
	BETA 200 S, BETA 200 M2, BETA 200 Ophthalmoscope	<b>X-01.88.069</b>	<b>X-02.88.070</b>
	BETA 200, 200 VET and K180 Otoscopes, Finoff, LAMBDA 100 Retinometer, Diascleral Transilluminator, Multi-purpose Illuminator, Slit Illumination Head	<b>X-01.88.077</b>	<b>X-02.88.078</b>
	Illumination adaptor for F.O. instruments		<b>X-02.88.078</b>
	DERMAPHOT 2.2 V vacuum lamp	<b>X-00.88.082</b>	
	K180 Ophthalmoscope	<b>X-01.88.084</b>	<b>X-02.88.086</b>
	BETA 200 Streak Retinoscope, alpha+ Streak Retinoscope	<b>X-01.88.087</b>	<b>X-02.88.089</b>
	BETA 200 Spot Retinoscope alpha+ Spot Retinoscope	<b>X-01.88.088</b>	<b>X-02.88.090</b>
	ClipLight vacuum lamp	<b>X-01.88.094</b>	
	HSL 100, HSL 150 and alpha+ HSL 150 Hand-held Slit Lamp	<b>X-01.88.098</b>	<b>X-02.88.099</b>

XHL Xenon Halogen bulbs are delivered in packs of 6; each bulb is individually blister-packed.