

## Universal Speculum



- ∴ **All-metal construction.** Opens symmetrically from 2.5 to 15mm.
- ∴ **Rotate ring to open branches.** Stays put in any position.
- ∴ **For use with Operating Otoscope, BETA100 Diagnostic Otoscope and Fiber Optic Nasal Illuminator.**

Universal Speculum

**B-000.11.239**

Not recommended for use with K100 Diagnostic Oscopes.

## Tip-Adaptor



- ∴ **All-metal construction.** Adaptor for AllSpec disposable tips (p. 15) for use with diagnostic (BETA100, K100) and operating otoscope.

Tip-Adaptor

**B-000.11.306**

## Insufflation Bulb



For Diagnostic Oscopes.

- ∴ **For pneumatic testing.**

Insufflation bulb for BETA100 and K100 Oscopes

**B-000.11.240**