

Each HEINE bulb is specifically designed as an integral part of the HEINE instrument it powers, based upon the following parameters:

Colour Temperature

Based on advanced HEINE XHL® Xenon Halogen Technology, each bulb is designed to provide the ideal colour temperature for the diagnosis the instrument is intended for. This ensures a proper balance between illumination brightness and true tissue colour for the most accurate diagnosis.

Bulb Life

HEINE bulbs last up to 50 % longer, providing long life, optimum performance and cost effectiveness.

Safety

Advanced design and manufacturing procedures prevent dangerous bulb implosions.

Optical Performance










The performance of each HEINE optical instrument depends on the bulb that powers it. Incorporated focusing lenses on the glass bulb envelope, precision filament placement, and homogeneous light projection are major factors in ensuring the satisfactory performance of your HEINE instrument.

We can only guarantee the performance and safety of your instrument if the appropriate HEINE bulb powers it. A non-genuine bulb or incorrect bulb will impede the performance of your HEINE instrument.

Insist on genuine HEINE replacement bulbs and parts to ensure the performance of your HEINE precision instrument.

HEINE XHL® Xenon Halogen Bulbs

For optimum performance and cost effectiveness

| | Instrument | 2.5 V | 3.5 V |
|---|---|----------------------------------|-------------|
|  | Autofoc, Examination Lamp with automatic connector | X-01.88.032 | X-02.88.047 |
|  | alpha Illuminated Loupe, mini 2000 Dermatoscope | X-01.88.034 | |
|  | alpha+ Illuminated Loupe, alpha+ Dermatoscope | X-01.88.034 | X-02.88.053 |
|  | Illumination head for Spreadable Anal Speculum, Glaucotest, F.O. Laryngoscope handle, angled F.O. Laryngoscope handle, Amnioscope – only 3.5 V, alpha+ Finoff | X-01.88.035 | X-02.88.044 |
|  | | F.O. SP Laryngoscope lamp handle | X-01.88.035 |
|  | Diagnostic Oscopes, K 100 and BETA 100, Operating Oscope, Curved Laryngeal Mirror (for 3.5 V: X-02.88.062), Tongue blade holder, alpha+ F.O. Oscope, alpha+ Tongue blade holder | X-01.88.037 | X-02.88.049 |
|  | mini 2000 F.O. Oscope, alpha Oscope, mini Tongue-blade holder, alpha Tongue-blade holder | X-01.88.037 | |
|  | Miroflex Ophthalmoscope, Illuminated Loupe with automatic connector, HFR 2 Spot Retinoscope, DELTA 10 | X-01.88.038 | X-02.88.046 |
|  | Focalux with automatic connector | X-01.88.039 | X-02.88.091 |
|  | Proctoscope head for proximal illumination | X-01.88.040 | X-02.88.051 |