



# **Dry Bath**

- Designed for Multi-functional Dry Block Heater
- User Self-Compensation Function to Control Difference between Real Temperature and Indicated Temperature (±10.0°C)
- Various kind of Blocks are offered (Optional)
- Over Temperature & Over Current Protection and Sensor Error Detection
- Patented Jog-Shuttle Control System and Digital LCD Display with Back-Light Function
- Weigh-performance Dry Bath-Heating ..... P. 43, 44

# >> Dry Bath - Heating & Cooling



Articles

# "witeg" WiseTherm® HB High-performance Dry Bath, Heating, up to $150^{\circ}$ C, $\pm 0.2^{\circ}$ C with Digital PID Control, Acrylic Lid, Precise Temperature Control, Modular Blocks, Anodized Aluminum Block

- Designed for Multi-functional Dry Block Heater
- Patented Jog-Shuttle Control System and Digital LCD Display with Back-Light Function
- **CE Certified**
- \* Ideal for Enzyme Reactions, Inactivation of Sera, Incubation and Other Laboratory Procedures
- \* Precise Temperature Control :  $\pm$  0,1 °C at 37 °C,  $\pm$  0,2 °C at 90 °C,  $\pm$ 0,4 °C at 135 °C
- \* User Self-Compensation Function to Control Difference between Real Temperature and Indicated Temperature (±10,0 °C)
- \* Over Temperature & Over Current Protection and Sensor Error Detection
- \* Digital Timer Included: up to 99 hr 59min (Operation Time + Delay Time)
- \* Storage Function of Temperature and Timer setting values
- \* Alarm Function: Error Status and Timer-end
- \* Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- \* The Block Extraction Tongs and Included for Easy Handing

#### Specification

Models		Single Control System		Dual Control System	
		"HB-48-Set"	"HB-96-Set"	"HB-96D-Set"	
Controller		Single Digital PID Controller Dual Digital PID Controller			
Capacity (1,5/2ml Tube)		48 Holes × 1,5/2ml Tube	es×1,5/2ml Tube 96 Holes×1,5/2ml Tube		
Heater		150W×2(Cartridge Heater)	150W×4 (Cartridge Heater)		
Dimension	Block (mm)	90×120×h60	(90×120×h60)×2ea	(90×120×h60)×2ea	
	Unit (mm)	180×220×h155	240×255×h155	285×250×h155	
Heater		150W×2(Cartridge Heater)	150W×4 (Ca	4 (Cartridge Heater)	
Temp. Range, & Accuracy		Ambient $+5$ °C to $150$ °C, $\pm 0.1$ °C at $37$ °C, $\pm 0.2$ °C at $90$ °C, $\pm 0.4$ °C at $135$ °C			
Heat-Up Speed & Sensor		Approx. 5℃/min, PT100			
Safety Circuit		Over Temp. and Over Current Protector, Sensor Error Detector			
Timer with Alarm		99hr 59min, Alarm (Error Status and Timer - end)			
Display		Digital LCD with Back-Light Function			
Material	Internal & External	Stainless Steel (#304), and Powder Coated Steel			
	Block	Anodized Aluminum			
Included Standard Acce.		Block "BLC548" for 1,5/2ml Tubes, and Block Tong "BLE101"			
Packing Size		400×270×h200 mm			
Gross Weight		4 kg	6 kg		
Power Consumption		300 W	600 W		
Power plug type		Hybrid E+F			
Power Supply*		AC 230V, 1N~, 50/60Hz			

<sup>\*</sup> Other Specifications are available upon Customer's Request.

### Ordering Information

# (a) Single Control System, available 1 Block

**DH.WHB00149** Dry Bath, Heating, **230V**, "**HB-48-Set**", with 1,5 / 2ml Block **DH.WHB00148** , "**HB-48 Unit**", without Block

#### (b) Single Control System, available 2 Blocks

**DH.WHB00197** Dry Bath, Heating, **230V**, "**HB-96-Set**", with 1,5 / 2ml Block **DH.WHB00196** , "**HB-96 Unit**", without Block

#### (c) Dual Control System, available 2 Blocks

**DH.WHB00199** Dry Bath, Heating, **230V**, "**HB-96D-Set**", with 1,5 / 2ml Block **DH.WHB00198** , "**HB-96D Unit**", without Block



(a) Single Control System, Dry Bath, Heating Package-Set "HB-48-Set" with Block "BLC548" for 1,5/2ml Centrifuge Tube ×48 holes



(b) Single Control System, Dry Bath, Heating Package-Set "HB-96-Set" with 2×Block "BLC548" for 1,5/2ml Centrifuge Tube ×48 holes



(c) Dual Control System, Dry Bath, Heating Package-Set "HB-96D-Set" with 2×Block "BLC548" for 1,5/2ml Centrifuge Tube ×48 holes

Articles

Cat. No

**DH.WHB535008**\* (10)

Description

# < Continued from... "witeg" WiseTherm® HB High-performance Dry Bath-Heating >



(a) Single Control System, Dry Bath, Heating Package-Set "HB-48-Set" with Block "BLC548" for 1,5/2ml Centrifuge Tube ×48 holes

#### Optional Interchangeable Blocks

DH.WHB510596\* (1) Block, "BLC596", for 0,2ml PCR Tube, 96 holes DH.WHB510570\* (2) Block, "BLC570", for 0,5ml Centrifuge Tube × 63 holes DH.WHB511548\* , "BLC548", for 1,5/2ml Centrifuge Tube × 48 holes DH.WHB521248\* (4) Block, "BLT248", for  $\phi$  12mm Test Tube, 48 holes DH.WHB521335\* , "BLT335", for  $\phi$  13mm Test Tube, 30 holes (5)DH.WHB521535\* "BLT535", for  $\phi$  15mm Test Tube, 30 holes DH.WHB521628\* , "BLT628", for  $\phi$  16,5mm Test Tube, 24 holes (7)DH.WHB521828\* , "BLT828", for  $\phi$  18mm Test Tube, 24 holes DH.WHB531515\* (9) Block, "BLC515", for 15ml Conical Tube, 15 holes

DH.WHB500101 Spare Block Extraction Tong, "BLE101"

\*1×Block Extraction Tongs, BLE101 included BLE101

, "BLC008", for 50ml Conical Tube, 8 holes


If you order BLT series Block, please make sure your test tube size is exact same with your order one.

#### Validation (IQ, OQ) Service


DH.WHB30148 Validation Service(IQ, OQ), "DBV148", for HB-48
DH.WHB30196 , "DBV196", for HB-96
DH.WHB30198 , "DBV198", for HB-96D

\* You will find matching tubes in our general catalogue at following pages: 223-225, 263-265, 371, 372.

# The Details of 10 Optional Blocks



**(1) "BLC596"** 0,2ml PCR Tube × 96 holes



**(2) "BLC570"** 0,5ml Centrifuge Tube × 63 holes



(3) "BLC548"

1,5 / 2ml Centrifuge Tube × 48 holes Included in the Package-Set



(4) "BLT248" *Ф* 12mm Tube × 48 holes



**(5) "BLT335"** *Ф* 13mm Tube × 30 holes



(6) "BLT535" *Ф* 15mm Tube × 30 holes



(7) "BLT628" *Φ* 16,5mm Tube × 24 holes



(8) "BLT828" *Ф* 18mm Tube × 24 holes



(9) "BLC515"
15ml Conical Tube × 15 holes



(10) "BLC008"
50ml Conical Tube × 8 holes

Description

Articles

# "witeg" WiseTherm® HB-R High-performance Dry Bath-Heating & Cooling with Precise Temperature Control, Modular Blocks, Acrylic Lid, Anodized Aluminum Block, −5 °C ~95 °C, ±0,15 °C

■ Designed for Multi-Function : Enzyme Reactions, Inactivation of Sera, Incubation and Other Laboratory Procedures.

- Patented Jog-Shuttle Control System and Digital LCD Display with Back-Light Function
- **CE Certified**
- \* User Self-Compensation Function to Control Difference between Real Temperature and Indicated Temperature
- \* Storage Function of Temperature and Timer setting values
- \* Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

#### Specification

Model		Heat/Cool Package-Set "HB - R48-Set"			
Capacity (1,5/2ml Tube)		48 holes × 1,5/2ml Centri. Tubes			
Dimension : Unit & Block		$185 \times 230 \times h220$ mm, and $90 \times 120 \times h60$ mm			
Cooling		Two Peltier Modules			
Cool-Down Speed		2,3°C/min. between 95°C to 25°C , 0,5°C/min. between 25°C to $-5$ °C			
Heater and Sensor		150 W × 2 (Cartridge Heater), PT100			
Temp. Range, & Accuracy		$-5^{\circ}\text{C}$ to $95^{\circ}\text{C}$ , and $\pm 0,1^{\circ}\text{C}$ at $37^{\circ}\text{C}$ , $\pm 0,2^{\circ}\text{C}$ at $90^{\circ}\text{C}$			
Heat-Up Speed		Approx. 5°C/min			
Safety Device		Over Temp. and Over Current Protector, Sensor Error Detector			
Timer with Alarm		99hr 59min, Alarm (Error Status and Timer - end)			
Display		Digital LCD Display with Back-Light Function			
Controller		Digital PID Controller			
Material	Internal & External	Stainless Steel (#304), and Powder Coated Steel			
	Block	Anodized Aluminum			
Included Standard Acce.		Block "BLC548" for 1,5/2ml Tubes, and Block Tong "BLE101"			
Packing Size, Gross Weight		350 × 220 × h310mm, 9 kg			
Power Consumption		450 W			
Power plug type		Hybrid E+F			
Power Supply*		AC 230V, 1N~, 50/60Hz			
Char Creatilestians are cyclickia way Quetamenta Degrees					

<sup>\*</sup> Other Specifications are available upon Customer's Request.

#### Ordering Information

## Dry Bath-Heating & Cooling Package-Set, with 1,5/2ml Block

DH.WHB00249 Dry Bath, Heat/Cool, "HB-R48-Set", 230V, with 1,5/2ml Block

#### Dry Bath-Heating & Cooling Unit Only, without Block

DH.WHB00248 Dry Bath, Heat/ Cool, "HB-R48-Unit", 230V, without Block

#### Optional Interchangeable Blocks

DH.WHB510596\* (1) Block, "BLC596", for 0,2ml PCR Tube, 96 holes DH.WHB510570\* (2) Block, "BLC570", for 0,5ml Centrifuge Tube × 63 holes DH.WHB511548\* , "BLC548", for 1,5/2ml Centrifuge Tube × 48 holes DH.WHB521248\* (4) Block, "BLT248", for  $\phi$  12mm Test Tube, 48 holes DH.WHB521335\* , "BLT335", for  $\phi$  13mm Test Tube, 30 holes DH.WHB521535\* (6) , "BLT535", for  $\phi$  15mm Test Tube, 30 holes "BLT628", for  $\phi$  16,5mm Test Tube, 24 holes DH.WHB521628\* DH.WHB521828\* , "BLT828", for  $\phi$  18mm Test Tube, 24 holes DH.WHB531515\* (9) Block, "BLC515", for 15ml Conical Tube, 15 holes DH.WHB535008\* , "BLC008", for 50ml Conical Tube, 8 holes

DH.WHB500101 Spare Block Extraction Tong, "BLE101"

 $*1 \times Block$  Extraction Tongs, BLE101 included

BLE101

If you order BLT series Block, please make sure your test tube size is exact same with your order one.

Validation (IQ, OQ) Service

DH.WHB30248 Validation Service(IQ, OQ), "DBRV248", for HB-R48

\* You will find matching tubes in our general catalogue at following pages: 223-225, 263-265, 371, 372.



Heat/Cool Package-Set "HB - R48-Set" with Block "BLC548" for 1,5/2ml Centrifuge Tube ×48 holes

# Included Block in the Package-Set



Articles

Cat. No

Description

# **Continued from...** "witeg" WiseTherm®HB-R High-performance Dry Bath-Heating & Cooling>

# The Details of 10 Optional Blocks



(1) "BLC596" 0,2ml PCR Tube × 96 holes



(2) "BLC570" 0,5ml Centrifuge Tube × 63 holes



(3) "BLC548"

1,5 / 2ml Centrifuge Tube × 48 holes Included in the Package-Set



**(4) "BLT248"** φ 12mm Tube × 48 holes



**(5) "BLT335"** *Ф* 13mm Tube × 30 holes



**(6)** "BLT535" *Ф* 15mm Tube × 30 holes



(7) "BLT628" *Φ* 16,5mm Tube × 24 holes



**(8) "BLT828"** *Ф* 18mm Tube × 24 holes



(9) "BLC515"
15ml Conical Tube × 15 holes



(10) "BLC008"
50ml Conical Tube × 8 holes