## DAkkS calibration

Without sticker

## Application:

Measuring and test equipment is frequently labelled with a calibration logo containing details of a calibration group. Use the table below to find and select the calibration services for your respective needs. The prices listed only apply to initial calibration of new products. When ordering, please enter your test equipment followed by the calibration service: Example of order: 1x no. 31342 025 **AURN**<sup>\*</sup> micrometer 25-50 mm incl. calibration

Art. no.	30000 505
Calibration Group	AN
Calibration Group Description	Precision pointer micrometers
Area of application for calibration group	25 - 100 mm
Gross Weight	0.001 kg
Product Group	387



