

Chirurgisches Nahtmaterial *Surgical Sutures*

SERAFAST®

PGACL

POLYGLYKOLSÄURE-CAPROLACTON
POLYGLYCOLIC ACID-CAPROLACTONE



violett / monofil / resorbierbar
violet / monofilament / absorbable

















ungefärbt / monofil / resorbierbar
undyed / monofilament / absorbable

Resorptionsprofil: 50 % Reißkraft nach 8 - 13 Tagen
0 % Reißkraft nach 90 - 120 Tagen

*Absorption profile: 50 % tensile strength after 8 - 13 days
0 % tensile strength after 90 - 120 days*



**SERAG
WIESSNER**

Nadelbezeichnung Needle		Stärke USP Size USP	Stärke metric Size metric	Fadenlänge/Meter Length of suture/metre	Abgabemenge/Stück Units per box	Bestell-Nummer ungefärbt Code number undyed	Bestell-Nummer violett Code number violet
DS-12		5/0	1	1 x 0,45 m	24	40101312	-
DS-15		4/0	1,5	1 x 0,70 m	24	40153213	-
DS-18		5/0	1	1 x 0,45 m	24	40101314	-
		5/0	1	1 x 0,70 m	24	40103214	JO103214
		4/0	1,5	1 x 0,45 m	24	40151314	-
		4/0	1,5	1 x 0,70 m	24	40153214	-
		3/0	2	1 x 0,45 m	24	40201314	-
		3/0	2	1 x 0,70 m	24	40203214	-
DS-25		4/0	1,5	1 x 0,70 m	24	40153215	-
		3/0	2	1 x 0,70 m	24	40203215	JO203215
		2/0	3	1 x 0,70 m	24	40303215	JO303215
DS-30		4/0	1,5	1 x 0,90 m	24	40154416	-
DSS-10		6/0	0,7	1 x 0,45 m	24	40071315	-
DSS-13		6/0	0,7	1 x 0,45 m	24	4007131L	-
		5/0	1	1 x 0,45 m	24	4010131L	JO10131L
		4/0	1,5	1 x 0,45 m	24	4015131L	-
DSS-18		4/0	1,5	1 x 0,45 m	24	4015131M	-
		4/0	1,5	1 x 0,70 m	24	4015321M	-
		3/0	2	1 x 0,45 m	24	4020131M	-
DSS-24		4/0	1,5	1 x 0,45 m	24	4015131N	-
		3/0	2	1 x 0,45 m	24	4020131N	-
		3/0	2	1 x 0,70 m	24	4020321N	-
FRX-27 +FRX-27		3/0	2	1 x 0,70 m	12	4P203226	-
GR-20		4/0	1,5	1 x 0,70 m	24	40153229	-
GS-12		4/0	1,5	1 x 0,45 m	24	40151333	-
GS-60		3/0	2	1 x 0,70 m	12	40203236	-


Nadelbezeichnung Needle		Stärke USP Size USP	Stärke metric Size metric	Fadenlänge/Meter Length of suture /metre	Abgabemenge/Stück Units per box	Bestell-Nummer ungefärbt Code number undyed	Bestell-Nummer violett Code number violet
HR-17		5/0	1	1 x 0,45 m	24	40101376	-
		5/0	1	1 x 0,70 m	24	40103276	JO103276
		4/0	1,5	1 x 0,70 m	24	40153276	JO153276
		3/0	2	1 x 0,70 m	24	40203276	JO203276
HR-17 + HR-17		4/0	1,5	1 x 0,70 m	36	-	JP153276
HR-22		4/0	1,5	1 x 0,70 m	24	4015322C	JO15322C
		3/0	2	1 x 0,70 m	24	4020322C	JO20322C
		3/0	2	4 x 0,45 m	12	-	JO20142C
		3/0	2	4 x 0,70 m	12	-	JO20412C
		2/0	3	1 x 0,70 m	24	-	JO30322C
HR-26		4/0	1,5	1 x 0,70 m	24	40153241	JO153241
		3/0	2	1 x 0,70 m	24	40203241	JO203241
		2/0	3	1 x 0,70 m	24	40303241	JO303241
HR-36		2/0	3	1 x 0,70 m	24	40303243	-


SERAFAST®

NADEL-FADEN-KOMBINATIONEN
NEEDLE SUTURE COMBINATIONS





SERAG-WIESSNER KG
Zum Kugelfang 8 - 12
95119 Naila/Germany


 + 49 9282 937-0


 + 49 9282 937-9369

Export Department:

 + 49 9282 9 37-230

 + 49 9282 937-9785

 info@serag-wiessner.de

 www.serag-wiessner.de