

DAkkS calibration

Without sticker

Application:

Measuring and test equipment is frequently labelled with a calibration logo containing details of a calibration group. Use the table below to find and select the calibration services for your respective needs. The prices listed only apply to initial calibration of new products. When ordering, please enter your test equipment followed by the calibration service: Example of order: 1x no. 31342 025 **ATORN** micrometer 25-50 mm incl. calibration

Art. no.	30000 561
Calibration Group	R
Calibration Group Description	Protractors
Gross Weight	0.001 kg
Product Group	387

